



The LikoLight is developed for those in need of a light mobile lift that folds away easily. LikoLight works much like any other mobile lift, but the lift is portable and can be conveniently stored. In fact, LikoLight can be quickly disassembled into three pieces to facilitate transport and storage.

Each detail on the LikoLight is carefully designed to provide maximum stability but minimum weight. The total weight is only 24 kg (53 lbs.) and no part weighs more than 10 kg (22 lbs). Practical carrier bags are available as an accessory.

Lifts with care

 **Liko**

## The versatile lightweight lift

Despite its size and weight, LikoLight is ideal for all common lifting situations; for example, lifts from bed to wheelchair, to and from the toilet and bathtub, or for lifting someone who has fallen to the floor.

LikoLight has electric raising and lowering and mechanical base-width adjustment. A patented easy-to-use mechanical emergency lowering is standard.

LikoLight Soft Bags are practical carrier bags made of a padded, durable and soft lightweight material. Perfect for travels.

### Important features:

- Folds away easily.
- Easily disassembled into three pieces.  
No part weighs more than 10 kg (22 lbs.).
- Made of aluminium; the right choice for environment, quality and design. The choice of material also makes the lift easy to roll and easy to keep clean.
- Mechanical emergency lowering.
- Stable construction.
- Ideal for traveling.
- Practical tabletop charger available as an accessory.



Accessory: LikoLight Soft bags



Easy to bring on the trip.



Manages all common lifting situations, including lifts from the floor.

### Technical data

Lifting capacity:	140 kg (308 lbs., 22 st.)
Lifting speed:	3.8 cm/s (1.5 inch/s) without load
Batteries:	2 x 12 V 2.9 Ah, valve regulated, sealed, lead accumulator (gel-type batteries)
Battery charger:	110-230 VAC 50-60 Hz, max. 400 mA
Material:	Aluminium
Weight:	24 kg (53 lbs.) total weight; disassembled, no part weighs more than 10 kg (22 lbs.)
Product numbers:	2030001 LikoLight 2037020 LikoLight Soft bags
Patent pending.	



LikoLight is tested by an accredited testing institute and complies with the requirements of directives for medical technical Class 1 products (MDD 93/42/EEC).

LikoLight complies with the requirements according to EN ISO 10535, IEC 60601-1, IEC 60601-1-2, UL-2601-1 and CAN/CSA C22.2 No.606.1.

Made in Sweden.

Liko is quality certified according to ISO 9001 and environmental certified according to ISO 14001.



# Liko

www.liko.com

Manufacturer:  
Liko AB  
SE-975 92 Luleå  
Sweden  
info@liko.se

Liko Distributor in UK:  
Liko (UK) Ltd • Brunel Centre  
Stroudwater Business Park, Stonehouse,  
Gloucester, GL10 3RU  
Tel: +44 (0)1453 827272  
Fax: +44 (0)1453 828844  
info@liko.co.uk

National Distributor in  
USA and Canada:  
Liko Inc • 842 Upper Union Street, Suite 4  
Franklin, MA 02038, USA  
Tel: 508-553-3993, 1-888-545 6671  
Fax: 508-528-6642  
info@likoinc.com