



Likorall Transfer Motor is intended for use with all Likorall lift motor models that are equipped with an outlet for connecting the transfer motor. This motor enables users to transfer themselves with the help of Liko along the rail. In such cases, the user should always have the use of a personal alarm to summon assistance in the event of an emergency. Another situation in which a Transfer motor is useful is when the lift must be moved over obstacles or limited-access areas.



Made in Sweden.

Liko is quality certified according to ISO 9001 and environmental certified according to ISO 14001.



Lifts with care

 Liko

### Technical data

Transfer speed: 15 cm/sec (5,9 in./sec)  
 Electronics: 24V 5A. Uses Liko's battery  
 Dimensions: 144x75x225 mm. WxHxL

Prod. No. 3126043

Reduces maximum lifting height by 50 mm.

[www.liko.com](http://www.liko.com)

Manufactured by:  
 Liko AB  
 Nedre vägen 100  
 SE-975 92 Luleå, Sweden  
[info@liko.se](mailto:info@liko.se)

National Distributor in USA and Canada:  
 Liko Inc  
 842 Upper Union Street, Suite 4  
 Franklin, MA 02038, USA  
 Phone 508-553-3993, 1-888-545-6671  
 Fax 508-528-6642  
[info@likoinc.com](mailto:info@likoinc.com)

Liko Distributor in UK:  
 Liko (UK) Ltd • Brunel Centre  
 Stroudwater Business Park, Stonehouse,  
 Gloucester, GL10 3RU  
 Tel: +44 (0)1453 827272  
 Fax: +44 (0)1453 828844  
[info@liko.co.uk](mailto:info@liko.co.uk)