

# Sabina II EE



# Sabina II EE

A flexible sit-to-stand lift intended for users who are able to actively participate in the standing exercise.



## Slingbars

For Sabina II you can choose between Sabina 350-slingbar and Comfort-slingbar.



## Seatstrap slingbar

Sabina II Seatstrap slingbar is mounted on the front of the slingbar and fits both the 350- and the Comfort-slingbar.



## Footrest

Sabina II has a new footrest prepared for heel-support and designed to prevent slipping.

Maximum load: 200 kg, 440 lbs

Made in Sweden

Liko is quality certified according to ISO 9002 and environment certified according to ISO 14001.

## Prod.no

2020003	Sabina II EE
2027002	Sabina II 350 slingbar
2027003	Sabina II Comfort slingbar
2027007	Sabina II Seatstrap slingbar
2027101	Sabina II Sidesupport

 **Liko**

[www.liko.com](http://www.liko.com)

Manufacturer

**Liko AB**

Nedre vägen 100 • SE-975 92 Luleå, Sweden

Tel +46 (0)920-24 47 00 • Fax +46 (0)920-24 47 01

[info@liko.se](mailto:info@liko.se)