

VIKING



Lifts with care



Just let your needs decide.

Viking is the result of 25 years of experience of lifting and transferring people in care settings.

All our knowledge of varying lifting needs, ergonomic and functional requirements, technical possibilities and modern manufacturing technology is incorporated in Viking's many features.

Together with Liko's comprehensive range of lifting accessories, Viking is your best friend when it comes to solving different lifting problems. Liko's slings and accessories, available in many different models and versions, will help you in virtually every situation, from light, easy lifts to lifting and transfer of patients weighing up to 880 lbs. (400 kg). For more information, read our brochure "Lifting Accessories".

Aluminum was the natural choice of material for Viking, not only for aesthetic reasons, but mainly because it makes your lift easier to maneuver. Aluminum also has a number of environmental advantages.

There are four Viking models: from super-strong XL to the smallest model, Viking S. A common feature of all models is the extremely good balance between low overall weight and high lifting capacity.

Just let your needs decide your choice!

VIKING S

VIKING M

VIKING L

VIKING XL

VIKING XL

For the most demanding lifts

Lifting heavy patients entails great responsibility. Although Viking XL has been specially developed for very heavy patients, it is also suitable for average-sized people. Despite its size, Viking XL weighs only 95 lbs. (43 kg) and is very easy to handle.

A lifting capacity of 660 lbs. (300 kg) makes Viking XL ideal for all heavy lifting situations, including lifts from the floor. Equipped with Armrests (accessory), is XL ready for walking training and can be an invaluable aid to rehabilitation, even for extremely heavy patients.

Heavy lifts are easy

Comfortable lifts up to 660 lbs. (300 kg) are no problem for Viking XL and Original Sling High Back.



Horizontal lifts in balance

With FlexoStretch and a Stretch leveller, even the heaviest patients can be lifted in a horizontal position.

Safe and secure gait-training

Armrests and Liko Lift Pants make Viking perfect for gait-training situations.



Viking XL

Lifting capacity: 660 lbs. (300 kg) Weight: 95 lbs. (43 kg) Art. No.: 2040003

VIKING **L**

Strong and easy to maneuver

Viking L is incredibly easy to maneuver in relation to its high lifting capacity and is ideal for all the most common lifting needs as well as for heavier patients. Like all Viking models, its low weight makes L very easy to handle.

Viking L can be equipped with Armrests, making it perfect for walking training. The product range includes accessories that have been designed specifically for walking training during rehabilitation.



All the way from the floor

With Original Sling High Back and Universal Twinbar 670, which provides extra space for the shoulders and body, heavy lifts from the floor are easily performed.



Down to the smallest detail

Large, easily maneuvered wheels are standard on Viking.

Perfect lifts

Safe and comfortable everyday-lifts with Liko Universal Sling.



Viking L

Lifting capacity: 550 lbs. (250 kg) Weight: 88 lbs. (40 kg) Art. No.: 2040004

VIKING **M**

The all-round lift

We are convinced that Viking M will become your new favourite lift. After 25 years of market-leading product development, we know how to design the optimal mobile lift.

Viking M has all of the best features and is ideal for the most common lifting situations.

Despite its low overall weight, Viking M can lift up to 450 lbs. (205 kg). This easy-to-manuevere lift even enables access to bathrooms where space is very limited.

Together with Liko's functional slings and accessories, Viking M is a perfect all-round lift for patients and caregivers with wide-ranging needs and high demands. Equipped with Armrests, Viking M is an invaluable aid to walking training and rehabilitation.



Easy gait-training

Liko MasterVest with front closure is easy to apply. With Armrests, Viking is ready for gait training.



In a horizontal position

Liko has many great accessories for horizontal lifts. Pictured is OctoStretch equipped with a Stretch leveller.

With extra care

The Cross-bar provides even more comfort for patients that are extra sensitive.



Viking M

Lifting capacity: 450 lbs. (205 kg) Weight: 77 lbs. (35 kg) Art. No.: 2040005

VIKING **S**

Small but strong

The smallest member of the Viking family is Viking S, a small, easily maneuvered mobile lift that performs superbly in relation to its diminutive size and low overall weight – perfect in settings where space is limited.

Viking S is appropriate for all common lifting situations; e.g., transfers between bed and wheelchair, to and from the toilet and bathtub or for lifting someone who has fallen to the floor.

Viking S weighs only 55 lbs. (25 kg) but can lift patients up to 350 lbs. (160 kg) – an unbeatable feat for a lift this size.

Lifts with Comfort

Viking S with Liko Original Sling, mod. 10, is a good combination for smooth lifts.



All the way from the floor

Original Sling High Back, mod. 20, provides a comfortable support during the lift.



Easy widening

Viking S has an easy-to-use mechanical widening of the base. Perfect when lifting from a wheelchair or the toilet.

Viking S

Lifting capacity: 350 lbs. (160 kg) Weight: 55 lbs. (25 kg) Art. No.: 2040006

Aluminum – the green metal



It wasn't simply by chance that Aluminum was the material chosen for the Viking series. Arguably, Aluminum production is energy-demanding, but the metal, once refined, will never rust. Low material weight makes the lifts easy to handle, with minimal load on the back and shoulders. And there are environmental benefits as well. Aluminum can be reused again and again, and recycling requires relatively little energy.

Lifting a person requires both technically good assistive devices and insight into medical implications. Liko is a specialist in the lifting and transfer of physically disabled persons. For more than 25 years, we have been successfully developing and manufacturing functional assistive devices that are safe and easy to use. Much of our success is based on a unique combination of engineering expertise, medical knowledge and our own experience of work in the care sector.

Liko is a Swedish company and the Nordic market leader. We are represented by subsidiaries and dealers in more than 25 countries worldwide.

Liko is accredited according to the ISO 9001 and ISO 14001 quality and environmental standards systems. Liko's lifts comply with the requirements of ISO 10535.



www.liko.com

Manufacturer:

Liko AB
SE-975 92 Luleå
Sweden
info@liko.se

National Distributor in USA and Canada:

Liko Inc • 842 Upper Union Street, Suite 4
Franklin, MA 02038, USA
Tel: 508-553-3993, 1-888-545 6671
Fax: 508-528-6642
info@likoinc.com

Liko Distributor in UK:

Liko (UK) Ltd • Brunel Centre
Stroudwater Business Park, Stonehouse,
Gloucester, GL10 3RU
Tel: +44 (0)1453 827272
Fax: +44 (0)1453 828844
info@liko.co.uk

www.liko.com