

PVI

MULTIFOLD RAMP

NEW AND IMPROVED
(Patent Pending)



This innovative multifold wheelchair ramp first introduced by PVI ten years ago now has several new features. Some of the improvements include a new design of the center joint which not only eliminates the pinch points of earlier models, it also adds strength to the highest stress area of the ramp. The ramp folds up and carries like a suitcase.

This new model is three inches shorter than the original model when folded making it easier to store and carry. Other improvements include auto locking of the

last two folds of the ramp so there is no need for a clip to hold the sides together when carrying. Just fold and go.

These ramps are all aluminum construction with a high traction tape on the surface. The multifold ramp is available in 5', 6', 7', 8', and 10' lengths. The ramp unfolds to a 30" wide platform that is ideal for wheelchairs and scooters with three or four wheels. The weight capacity is 600 pounds. You can expect quick delivery by UPS of these stock ramps. Custom widths are available.

"Multifold Ramp" - Specifications

ITEM NUMBER	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	CAPACITY
WCR530	5 feet	30 inches	32½" x 16" x 8½"	27 lbs.	600 lbs.
WCR630	6 feet	30 inches	38½" x 16" x 8½"	32 lbs.	600 lbs.
WCR730	7 feet	30 inches	44½" x 16" x 8½"	38 lbs.	600 lbs.
WCR830	8 feet	30 inches	50½" x 16" x 9"	45 lbs.	600 lbs.
WC1030	10 feet	30 inches	62½" x 16" x 9"	56 lbs.	600 lbs.