



SINGLEFOLD RAMP



The Singlefold Ramp, designed by PVI, is easy to handle and set up. The ramp folds in the middle and has a convenient carrying handle. The ramp unfolds to a 30 inch wide platform that is ideal for wheelchairs and scooters with three or four wheels. These

single fold ramps are the lightest weight on the market. Other features include all aluminum construction and high traction tape on the surface. The weight capacity is 600 pounds. You can expect quick delivery by UPS of this stock item. Custom widths are available.

"Singlefold Ramp" - Specifications

ITEM NUMBER	LENGTH	WIDTH	FOLDED SIZE	WEIGHT	CAPACITY
SFW230	2 feet	30 inches	24" x 16" x 4"	10 lbs.	600 lbs.
SFW330	3 feet	30 inches	36" x 16" x 4"	14 lbs.	600 lbs.
SFW430	4 feet	30 inches	48" x 16" x 4"	18 lbs.	600 lbs.
SFW530	5 feet	30 inches	60" x 16" x 4"	22 lbs.	600 lbs.