

# PVI

# SOLID

The PVI ramp is an excellent choice when a more permanent ramp is desired. No assembly required. Its one-piece design makes it ideal for wheelchairs and scooters with three or four wheels. This ramp is constructed of aluminum and has a slip resistant surface. The PVI solid ramp is available in three length and two widths. Custom widths are also available. The load capacity is 600 pounds.



### “Solid” - Specifications

ITEM NUMBER	LENGTH	WIDTH	WEIGHT	CAPACITY
SL330	3 feet	30 inches	13 lbs.	600 lbs.
SL336	3 feet	36 inches	19 lbs.	600 lbs.
SL430	4 feet	30 inches	18 lbs.	600 lbs.
SL436	4 feet	36 inches	24 lbs.	600 lbs.
SL530	5 feet	30 inches	23 lbs.	600 lbs.
SL536	5 feet	36 inches	29 lbs.	600 lbs.