

# Reliant™

## Reliability in a Compact Design

- Compact Design
- Automatic Inboard - Outboard Rollstops
- Quiet, Smooth Operation
- 5-Year Limited Warranty
- Easy to Maintain
- Affordably Priced



**T**he automatic Reliant wheelchair lift from Ricon offers you dependable accessibility while its trim, compact design maximizes valuable interior space.

**The Reliant also features a strong, durable steel frame built to give you years of dependable service.** Its patented Sto-Loc™ design assures a quiet interior ride by securely locking the platform in the stowed position.

So give yourself the lift that features our very best—the space saving, affordable, reliable Reliant wheelchair lift from Ricon.



U.S. Patents: Nos. 5,308,215; 5,445,488; 5,605,431 and 5,944,473; Australian Patent Nos. 661127; 687066; 705929; Canadian Patent No. 2,168,761; other U.S. and foreign patents pending. ©Copyright 2001, Ricon Corp. Printed in U.S.A.

**RICON®**

*Innovation in Mobility*

# Reliant™

## Standard Features

- 31" (787 mm) Lifting Height
- Convenient Handrails for User Comfort and Security
- Non-Skid Platform
- 5-Year Limited Warranty
- Welded All-Steel Frame
- Automatic Inboard and Outboard Rollstops
- Powerful Hydraulic Pump
- Built-In Manual Backup System
- Simple Electrical System
- Permanent, No-Lube Bearings at All Major Pivot Points




*Your vehicle can be equipped with a Ricon automatic door operator for complete independence and convenience.*

### Specifications

Power	Electro - Hydraulic Motor
Pump Rating	12 Volt DC, 65 Amp Avg. / Cycle, 1250 psi / 8.6 MPa
Hydraulic Cylinder(s)	(2) 1.5" (38mm), Power Up / Power Down
Rated Load Capacity	800 Lbs. (363 Kg)
Manual Backup - Up	Manual Hand Pump
Manual Backup - Down	Manual Pressure Release Valve
Weight	317 Lbs. (144 Kg)

### Your Choice of Options Includes:

- A. Internal Swing-Door Operator
- B. Slide-Door Operator
- C. Remote Control



A



B



C

**RICON®**

*Innovation in Mobility*